

For Research Use Only

CHD6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83122-5-PBS

Basic Information

Catalog Number: 83122-5-PBS	GenBank Accession Number: NM_032221	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 84181	CloneNo.: 230397C3
Source: Rabbit	UNIPROT ID: Q8TD26	
Isotype: IgG	Full Name: chromodomain helicase DNA binding protein 6	
	Calculated MW: 305 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00298-3, CHD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83122-6-PBS. Detection antibody: 83122-5-PBS. Standard: SY01008. Range: 0.156-5 ng/mL