

For Research Use Only

# NLRP12 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83123-1-PBS

## Basic Information

<b>Catalog Number:</b> 83123-1-PBS	<b>GenBank Accession Number:</b> BC028069	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 91662	<b>CloneNo.:</b> 230404D9
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P59046	
<b>Isotype:</b> IgG	<b>Full Name:</b> NLR family, pyrin domain containing 12	
<b>Immunogen Catalog Number:</b> AG34281	<b>Calculated MW:</b> 1061 aa, 120 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00300-3, NLRP12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83123-5-PBS. Detection antibody: 83123-1-PBS. Standard: Ag34281. Range: 2.5-40 ng/mL.