

For Research Use Only

NLRP12 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83123-3-PBS

Basic Information

Catalog Number: 83123-3-PBS	GenBank Accession Number: BC028069	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 91662	CloneNo.: 230404B8
Source: Rabbit	UNIPROT ID: P59046	
Isotype: IgG	Full Name: NLR family, pyrin domain containing 12	
Immunogen Catalog Number: AG34281	Calculated MW: 1061 aa, 120 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00300-1, NLRP12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83123-3-PBS. Detection antibody: 83123-2-PBS. Standard: Ag34281. Range: 2.5-40 ng/mL.



Cytometric bead array standard curve of MP00300-2, NLRP12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83123-3-PBS. Detection antibody: 83123-4-PBS. Standard: Ag34281. Range: 2.5-40 ng/mL.