

For Research Use Only

CDX2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83125-4-PBS

Basic Information

Catalog Number:

83125-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001265

GeneID (NCBI):

1045

UNIPROT ID:

Q99626

Full Name:

caudal type homeobox 2

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

230289G12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00304-1, CDX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83125-1-PBS. Detection antibody: 83125-4-PBS. Standard: SY01336. Range: 0.313-10 ng/mL



Cytometric bead array standard curve of MP00304-2, CDX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83125-3-PBS. Detection antibody: 83125-4-PBS. Standard: SY01336. Range: 0.313-10 ng/mL