

For Research Use Only

FOXD1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83126-3-PBS

Basic Information

Catalog Number: 83126-3-PBS	GenBank Accession Number: NM_004472.2	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2297	CloneNo.: 230377F7
Source: Rabbit	UNIPROT ID: Q16676	
Isotype: IgG	Full Name: forkhead box D1	
Immunogen Catalog Number: AG29344	Calculated MW: 46 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00301-1, FOXD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83126-2-PBS. Detection antibody: 83126-3-PBS. Standard: Ag29344. Range: 0.625-20 ng/mL