

For Research Use Only

GRIA2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83128-1-PBS

Basic Information

Catalog Number: 83128-1-PBS	GenBank Accession Number: BC010574	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2891	CloneNo.: 230403D1
Source: Rabbit	UNIPROT ID: P42262	
Isotype: IgG	Full Name: glutamate receptor, ionotropic, AMPA 2	
Immunogen Catalog Number: AG2662	Calculated MW: 883 aa, 99 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00143-1, GRIA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83128-2-PBS. Detection antibody: 83128-1-PBS. Standard: Ag2662. Range: 1.25-20 ng/mL