

For Research Use Only

HNRPDL Recombinant antibody, PBS Only (Detector)



Catalog Number: 83129-3-PBS

Basic Information

Catalog Number: 83129-3-PBS	GenBank Accession Number: BC007392	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9987	CloneNo.: 230415B3
Source: Rabbit	UNIPROT ID: O14979	
Isotype: IgG	Full Name: heterogeneous nuclear ribonucleoprotein D-like	
Immunogen Catalog Number: AG1046	Calculated MW: 46 aa, 4 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00151-1, HNRPDL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83129-4-PBS. Detection antibody: 83129-3-PBS. Standard: Ag1046. Range: 1.25-20 ng/mL.