

For Research Use Only

MTCO2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83132-4-PBS

Basic Information

Catalog Number: 83132-4-PBS	GenBank Accession Number: YP_003024029	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4513	CloneNo.: 230448C11
Source: Rabbit	Full Name: cytochrome c oxidase II	
Isotype: IgG	Calculated MW: 26 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00254-3, MTCO2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83132-4-PBS. Detection antibody: 83132-6-PBS. Standard: SY01099. Range: 0.313-20 ng/mL