

For Research Use Only

ORM2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83135-1-PBS

Basic Information

Catalog Number: 83135-1-PBS	GenBank Accession Number: BC015964	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5005	CloneNo.: 230477B2
Source: Rabbit	UNIPROT ID: P19652	
Isotype: IgG	Full Name: orosomuroid 2	
Immunogen Catalog Number: AG1667	Calculated MW: 24 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00283-2, ORM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83135-1-PBS. Detection antibody: 83135-3-PBS. Standard: Ag1667. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00283-3, ORM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83135-1-PBS. Detection antibody: 83135-2-PBS. Standard: Ag1667. Range: 1.25-20 ng/mL