

For Research Use Only

# RNF31 / HOIP Recombinant antibody



Catalog Number: 83137-4-RR

## Basic Information

Catalog Number:

83137-4-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34505

GenBank Accession Number:

BC017376

GeneID (NCBI):

55072

UNIPROT ID:

Q96EPO

Full Name:

ring finger protein 31

Calculated MW:

1072 aa, 120 kDa

Observed MW:

100 kDa

Purification Method:

Protein A purification

CloneNo.:

230446E9

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HEK-293 cells, HeLa cells, Jurkat cells, MCF-7 cells, HepG2 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83137-4-RR (RNF31 / HOIP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.