

For Research Use Only

# CTNND2 Recombinant antibody, PBS Only



Catalog Number: 83138-5-PBS

## Basic Information

Catalog Number:  
83138-5-PBS

Size:  
1mg/ml

Source:  
Rabbit

Isotype:  
IgG

GenBank Accession Number:  
NM\_001332

GeneID (NCBI):  
1501

UNIPROT ID:  
Q9UQB3

Full Name:  
catenin (cadherin-associated protein),  
delta 2 (neural plakophilin-related  
arm-repeat protein)

Calculated MW:  
133 kDa

Purification Method:  
Protein A purification

CloneNo.:  
230396F4

## Applications

Tested Applications:  
WB, Indirect ELISA

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



rat brain tissue were subjected to SDS PAGE followed by western blot with 83138-5-RR (CTNND2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83138-5-PBS in a different storage buffer formulation.