

For Research Use Only

TGS1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83140-4-PBS

Basic Information

Catalog Number: 83140-4-PBS	GenBank Accession Number: BC011999	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 96764	CloneNo.: 230427G7
Source: Rabbit	UNIPROT ID: Q96R50	
Isotype: IgG	Full Name: trimethylguanosine synthase homolog (<i>S. cerevisiae</i>)	
Immunogen Catalog Number: AG3683	Calculated MW: 852 aa, 90 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00299-1, TGS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83140-2-PBS. Detection antibody: 83140-4-PBS. Standard: Ag3683. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00299-2, TGS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83140-3-PBS. Detection antibody: 83140-4-PBS. Standard: Ag3683. Range: 0.625-20 ng/mL