

For Research Use Only

MDM1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83141-1-PBS

Basic Information

Catalog Number: 83141-1-PBS	GenBank Accession Number: BC028355	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 56890	CloneNo.: 230425A4
Source: Rabbit	UNIPROT ID: Q8TC05	
Isotype: IgG	Full Name: Mdm1 nuclear protein homolog (mouse)	
Immunogen Catalog Number: AG10541	Calculated MW: 81 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00166-2, MDM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83141-1-PBS. Detection antibody: 83141-2-PBS. Standard: Ag10541. Range: 0.625-20 ng/mL