

For Research Use Only

P2RY1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83143-1-PBS

Basic Information

Catalog Number: 83143-1-PBS	GenBank Accession Number: BC074785	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5028	CloneNo.: 230473C3
Source: Rabbit	UNIPROT ID: P47900	
Isotype: IgG	Full Name: purinergic receptor P2Y, G-protein coupled, 1	
Immunogen Catalog Number: AG12992	Calculated MW: 373 aa, 42 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00175-1, P2RY1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83143-2-PBS. Detection antibody: 83143-1-PBS. Standard: Ag12992. Range: 0.625-20 ng/mL