

For Research Use Only

COQ4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83144-1-PBS

Basic Information

Catalog Number: 83144-1-PBS	GenBank Accession Number: BC011895	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 51117	CloneNo.: 230482E8
Source: Rabbit	UNIPROT ID: Q9Y3A0	
Isotype: IgG	Full Name: coenzyme Q4 homolog (S. cerevisiae)	
Immunogen Catalog Number: AG10089	Calculated MW: 265 aa, 30 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00179-1, COQ4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83144-2-PBS. Detection antibody: 83144-1-PBS. Standard: Ag10089. Range: 0.625-20 ng/mL