

For Research Use Only

GBX2 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83145-4-PBS

Basic Information

Catalog Number: 83145-4-PBS	GenBank Accession Number: BC137448	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 2637	CloneNo.: 230468G8
Source: Rabbit	UNIPROT ID: P52951	
Isotype: IgG	Full Name: gastrulation brain homeobox 2	
Immunogen Catalog Number: AG16288	Calculated MW: 348 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00200-1, GBX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83145-2-PBS. Detection antibody: 83145-4-PBS. Standard: Ag16288. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00200-2, GBX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83145-3-PBS. Detection antibody: 83145-4-PBS. Standard: Ag16288. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00200-3, GBX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83145-1-PBS. Detection antibody: 83145-4-PBS. Standard: Ag16288. Range: 0.625-20 ng/mL