

For Research Use Only

GRP94 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83147-1-PBS

Basic Information

Catalog Number:

83147-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1439

GenBank Accession Number:

BC009195

GeneID (NCBI):

7184

UNIPROT ID:

P14625

Full Name:

heat shock protein 90kDa beta (Grp94), member 1

Calculated MW:

96 kDa

Purification Method:

Protein A purification

CloneNo.:

230502D10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00194-3, GRP94 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83147-2-PBS. Detection antibody: 83147-1-PBS. Standard: Ag1439. Range: 0.625-20 ng/mL