

For Research Use Only

CD276 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83152-4-PBS

Basic Information

Catalog Number:

83152-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001024736

GeneID (NCBI):

80381

UNIPROT ID:

Q5ZPR3

Full Name:

CD276 molecule

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

230305F11

Applications

Tested Applications:

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00248-1, Human B7-H3/CD276 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83152-4-PBS. Detection antibody: 83152-1-PBS. Standard: Eg0075. Range: 0.625-20 ng/mL.