

For Research Use Only

FGFR4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83153-1-PBS

Basic Information

| | | |
|---------------------------------------|--|---|
| Catalog Number: 83153-1-PBS | GenBank Accession Number: BC011847 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 2264 | CloneNo.: 230364C11 |
| Source: Rabbit | UNIPROT ID: P22455 | |
| Isotype: IgG | Full Name: fibroblast growth factor receptor 4 | |
| | Calculated MW: 88 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00257-1, FGFR4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83153-4-PBS. Detection antibody: 83153-1-PBS. Standard: Eg0221. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00257-2, FGFR4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83153-2-PBS. Detection antibody: 83153-1-PBS. Standard: Eg0221. Range: 0.625-20 ng/mL