

For Research Use Only

# FGFR4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83153-3-PBS

## Basic Information

Catalog Number:

83153-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC011847

GeneID (NCBI):

2264

UNIPROT ID:

P22455

Full Name:

fibroblast growth factor receptor 4

Calculated MW:

88 kDa

Purification Method:

Protein A purification

CloneNo.:

230364D9

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00257-3, FGFR4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83153-2-PBS. Detection antibody: 83153-3-PBS. Standard: Eg0221. Range: 0.625-20 ng/mL