

For Research Use Only

HOXB6 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83155-1-PBS

Basic Information

Catalog Number: 83155-1-PBS	GenBank Accession Number: BC014651	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3216	CloneNo.: 230433C11
Source: Rabbit	UNIPROT ID: P17509	
Isotype: IgG	Full Name: homeobox B6	
Immunogen Catalog Number: AG14189	Calculated MW: 224 aa, 25 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00150-1, HOXB6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83155-1-PBS. Detection antibody: 83155-2-PBS. Standard: Ag14189. Range: 1.25-20 ng/mL



Cytometric bead array standard curve of MP00150-2, HOXB6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83155-1-PBS. Detection antibody: 83155-3-PBS. Standard: Ag14189. Range: 2.5-20 ng/mL