

For Research Use Only

# Mouse MCP-1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83157-1-PBS

## Basic Information

Catalog Number:

83157-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM-011333

GeneID (NCBI):

20296

Full Name:

chemokine (C-C motif) ligand 2

Purification Method:

Protein A purification

CloneNo.:

230483C2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00168-3, Mouse Ccl2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83157-4-PBS. Detection antibody: 83157-1-PBS. Standard: Eg0427. Range: 1.25-20 ng/mL