

For Research Use Only

PD-L2/B7-DC Recombinant antibody, PBS Only (Capture)

Catalog Number: 83158-2-PBS



Basic Information

Catalog Number:

83158-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074766

GeneID (NCBI):

80380

UNIPROT ID:

Q9BQ51

Full Name:

programmed cell death 1 ligand 2

Calculated MW:

273 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

230514H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00778-2, PD-L2/B7-DC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83158-2-PBS. Detection antibody: 83158-4-PBS. Standard: Eg0507. Range: 0.625-20 ng/mL.