

For Research Use Only

PD-L2/B7-DC Recombinant antibody, PBS Only (Detector)

Catalog Number: 83158-4-PBS



Basic Information

Catalog Number: 83158-4-PBS	GenBank Accession Number: BC074766	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 80380	CloneNo.: 230514H6
Source: Rabbit	UNIPROT ID: Q9BQ51	
Isotype: IgG	Full Name: programmed cell death 1 ligand 2	
	Calculated MW: 273 aa, 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00778-2, PD-L2/B7-DC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83158-2-PBS. Detection antibody: 83158-4-PBS. Standard: Eg0507. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00778-3, PD-L2/B7-DC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83158-3-PBS. Detection antibody: 83158-4-PBS. Standard: Eg0507. Range: 0.625-20 ng/mL.