

For Research Use Only

# TNFAIP6 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83163-2-PBS



## Basic Information

Catalog Number:

83163-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3885

GenBank Accession Number:

BC030205

GeneID (NCBI):

7130

UNIPROT ID:

P98066

Full Name:

tumor necrosis factor, alpha-induced protein 6

Calculated MW:

277 aa, 31 kDa

Purification Method:

Protein A purification

CloneNo.:

230438E11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00218-1, TNFAIP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83163-2-PBS. Detection antibody: 83163-3-PBS. Standard: Ag3885. Range: 0.625-20 ng/mL.