

For Research Use Only

Mouse Angiopoietin-like 3/ANGPTL3 Recombinant antibody, PBS Only



Catalog Number: 83165-3-PBS

Basic Information

Catalog Number: 83165-3-PBS	GenBank Accession Number: NM_013913.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 30924	CloneNo.: 230498F11
Source: Rabbit	UNIPROT ID: Q9R182	
Isotype: IgG	Full Name: angiopoietin-like 3	
	Calculated MW: 53KD	

Applications

Tested Applications:
WB, ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83165-3-RR (Angiopoietin-like 3/ANGPTL3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83165-3-PBS in a different storage buffer formulation.



A375 cells were subjected to SDS PAGE followed by western blot with 83165-3-RR (Angiopoietin-like 3/ANGPTL3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83165-3-PBS in a different storage buffer formulation.