

For Research Use Only

ROCK1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83166-2-PBS

Basic Information

Catalog Number: 83166-2-PBS	GenBank Accession Number: BC113114	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 6093	CloneNo.: 230423E7
Source: Rabbit	UNIPROT ID: Q13464	
Isotype: IgG	Full Name: Rho-associated, coiled-coil containing protein kinase 1	
Immunogen Catalog Number: AG16317	Calculated MW: 1354 aa, 158 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00229-1, ROCK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83166-4-PBS. Detection antibody: 83166-2-PBS. Standard: Ag16317. Range: 0.625-20 ng/mL