

For Research Use Only

TGFB2-Specific Recombinant antibody, PBS Only (Capture)



Catalog Number:83167-1-PBS

Basic Information

Catalog Number:

83167-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003238

GeneID (NCBI):

7042

UNIPROT ID:

P61812

Full Name:

transforming growth factor, beta 2

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230461C9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00259-3, TGFβ2-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83167-1-PBS. Detection antibody: 83167-2-PBS. Standard: SY00619. Range: 0.156-10 ng/mL.