

For Research Use Only

TGFB2-Specific Recombinant antibody, PBS Only



Catalog Number: 83167-5-PBS

Basic Information

Catalog Number:

83167-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003238

GeneID (NCBI):

7042

UNIPROT ID:

P61812

Full Name:

transforming growth factor, beta 2

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

230461F5

Applications

Tested Applications:

IF, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using TGF β 2-Specific antibody (83167-5-RR, Clone: 230461F5) at dilution of 1:250 and CoraLite $\text{\textcircled{R}}$ 488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1). This data was developed using the same antibody clone with 83167-5-PBS in a different storage buffer formulation.