

For Research Use Only

SAR1A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83168-4-PBS

Basic Information

Catalog Number: 83168-4-PBS	GenBank Accession Number: NM_020150	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 56681	CloneNo.: 230471C1
Source: Rabbit	UNIPROT ID: Q9NR31	
Isotype: IgG	Full Name: SAR1 homolog A (S. cerevisiae)	
	Calculated MW: 22 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00263-2, SAR1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83168-1-PBS. Detection antibody: 83168-4-PBS. Standard: SY01131. Range: 0.078-5 ng/mL



Cytometric bead array standard curve of MP00263-3, SAR1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83168-5-PBS. Detection antibody: 83168-4-PBS. Standard: SY01131. Range: 0.156-5 ng/mL