

For Research Use Only

SMAD4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83169-4-PBS

Basic Information

Catalog Number: 83169-4-PBS	GenBank Accession Number: BC002379	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4089	CloneNo.: 240026E5
Source: Rabbit	UNIPROT ID: Q13485	
Isotype: IgG	Full Name: SMAD family member 4	
Immunogen Catalog Number: AG0299	Calculated MW: 60 kDa	

Applications

Tested Applications:
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00267-1, Human SMAD4 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83169-4-PBS. Detection antibody: 83169-2-PBS. Standard: Ag0299. Range: 0.625-20 ng/mL



Standard curve of MP00267-2, Human SMAD4 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83169-4-PBS. Detection antibody: 83169-1-PBS. Standard: Ag0299. Range: 0.625-20 ng/mL