

For Research Use Only

# Histone H2B Recombinant antibody, PBS Only (Capture)



Catalog Number: 83171-1-PBS

## Basic Information

Catalog Number:

83171-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC005827

GeneID (NCBI):

8349

UNIPROT ID:

Q16778

Full Name:

histone cluster 2, H2be

Calculated MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

230539A10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00183-2, Histone H2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83171-1-PBS. Detection antibody: 83171-3-PBS. Standard: SY02145Ac. Range: 0.625-5 ng/mL