

For Research Use Only

# USP8 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83172-3-PBS

## Basic Information

Catalog Number: 83172-3-PBS	GenBank Accession Number: BC110590	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9101	CloneNo.: 240016F7
Source: Rabbit	UNIPROT ID: P40818	
Isotype: IgG	Full Name: ubiquitin specific peptidase 8	
Immunogen Catalog Number: AG27143	Calculated MW: 1118 aa, 128 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00188-2, USP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83172-2-PBS. Detection antibody: 83172-3-PBS. Standard: Ag27143. Range: 0.625-20 ng/mL