

For Research Use Only

# TRIM54 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83175-3-PBS

## Basic Information

Catalog Number: 83175-3-PBS	GenBank Accession Number: BC141807	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 57159	CloneNo.: 230391F2
Source: Rabbit	UNIPROT ID: Q9BYV2	
Isotype: IgG	Full Name: tripartite motif-containing 54	
Immunogen Catalog Number: AG15229	Calculated MW: 358 aa, 40 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00236-2, TRIM54 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83175-3-PBS. Detection antibody: 83175-4-PBS. Standard: Ag15229. Range: 0.625-20 ng/mL