

For Research Use Only

# Rat IL-6 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83179-1-PBS

## Basic Information

Catalog Number:

83179-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_012589.2

GeneID (NCBI):

24498

UNIPROT ID:

P20607

Full Name:

interleukin 6

Calculated MW:

24 kDa

Purification Method:

Protein A purification

CloneNo.:

240003B1

## Applications

Tested Applications:

Species Specificity:

Rat

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00264-1, Rat IL-6 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83179-1-PBS. Detection antibody: 83179-3-PBS. Standard: Eg0800. Range: 0.625-20 ng/mL



Standard curve of MP00264-3, Rat IL-6 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83179-2-PBS. Detection antibody: 83179-1-PBS. Standard: Eg0800. Range: 0.625-20 ng/mL