

For Research Use Only

IL11RA Recombinant antibody, PBS Only (Detector)



Catalog Number: 83180-2-PBS

Basic Information

Catalog Number: 83180-2-PBS	GenBank Accession Number: BC003110	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3590	CloneNo.: 240019F3
Source: Rabbit	UNIPROT ID: Q14626	
Isotype: IgG	Full Name: interleukin 11 receptor, alpha	
Immunogen Catalog Number: AG28360	Calculated MW: 45 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00268-1, IL11RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83180-4-PBS. Detection antibody: 83180-2-PBS. Standard: Ag28360. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00268-2, IL11RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83180-3-PBS. Detection antibody: 83180-2-PBS. Standard: Ag28360. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00268-3, IL11RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83180-1-PBS. Detection antibody: 83180-2-PBS. Standard: Ag28360. Range: 0.625-20 ng/mL