

For Research Use Only

# BBS13 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83181-3-PBS

## Basic Information

<b>Catalog Number:</b> 83181-3-PBS	<b>GenBank Accession Number:</b> BC010061	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 54903	<b>CloneNo.:</b> 240036E2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9NXB0	
<b>Isotype:</b> IgG	<b>Full Name:</b> Meckel syndrome, type 1	
<b>Immunogen Catalog Number:</b> AG9177	<b>Calculated MW:</b> 559 aa, 65 kDa	

## Applications

**Tested Applications:**  
**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00270-2, Human BBS13 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83181-2-PBS. Detection antibody: 83181-3-PBS. Standard: Ag9177. Range: 0.625-20 ng/mL



Standard curve of MP00270-3, Human BBS13 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83181-4-PBS. Detection antibody: 83181-3-PBS. Standard: Ag9177. Range: 0.625-20 ng/mL