

For Research Use Only

BBS13 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83181-4-PBS

Basic Information

Catalog Number: 83181-4-PBS	GenBank Accession Number: BC010061	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 54903	CloneNo.: 240036G2
Source: Rabbit	UNIPROT ID: Q9NXB0	
Isotype: IgG	Full Name: Meckel syndrome, type 1	
Immunogen Catalog Number: AG9177	Calculated MW: 559 aa, 65 kDa	

Applications

Tested Applications:
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00270-3, Human BBS13 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83181-4-PBS. Detection antibody: 83181-3-PBS. Standard: Ag9177. Range: 0.625-20 ng/mL.