

For Research Use Only

SLC38A9 Recombinant antibody, PBS Only



Catalog Number: 83182-5-PBS

Basic Information

Catalog Number: 83182-5-PBS	GenBank Accession Number: BC066891	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 153129	CloneNo.: 240039E7
Source: Rabbit	UNIPROT ID: Q8NBW4	
Isotype: IgG	Full Name: solute carrier family 38, member 9	
Immunogen Catalog Number: AG15607	Calculated MW: 561 aa, 64 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse pancreas tissue was subjected to SDS PAGE followed by western blot with 83182-5-RR (SLC38A9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83182-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83182-5-RR against Human SLC38A9 were performed. The affinity constant is 0.715 nM.