

For Research Use Only

Hemoglobin Alpha Recombinant antibody, PBS Only (Capture)



Catalog Number: 83185-2-PBS

Basic Information

Catalog Number: 83185-2-PBS	GenBank Accession Number: BC050661	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3039	CloneNo.: 240049G8
Source: Rabbit	UNIPROT ID: P69905	
Isotype: IgG	Full Name: hemoglobin, alpha 1	
Immunogen Catalog Number: AG6040	Calculated MW: 15 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00239-2, HBA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83185-2-PBS. Detection antibody: 83185-3-PBS. Standard: Ag6040. Range: 3.125-100 ng/mL.