

For Research Use Only

SMC6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83187-1-PBS

Basic Information

Catalog Number: 83187-1-PBS	GenBank Accession Number: BC039828	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 79677	CloneNo.: 240015C2
Source: Rabbit	UNIPROT ID: Q965B8	
Isotype: IgG	Full Name: structural maintenance of chromosomes 6	
Immunogen Catalog Number: AG6898	Calculated MW: 126 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00255-1, SMC6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83187-2-PBS. Detection antibody: 83187-1-PBS. Standard: Ag6898. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00255-2, SMC6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83187-3-PBS. Detection antibody: 83187-1-PBS. Standard: Ag6898. Range: 1.25-20 ng/mL