

For Research Use Only

MMP13 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83188-1-PBS

Basic Information

Catalog Number: 83188-1-PBS	GenBank Accession Number: BC074808	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4322	CloneNo.: 240060A6
Source: Rabbit	UNIPROT ID: P45452	
Isotype: IgG	Full Name: matrix metalloproteinase 13	
Immunogen Catalog Number: AG34892	Calculated MW: 471 aa, 54 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00266-2, Human MMP13 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83188-1-PBS. Detection antibody: 83188-2-PBS. Standard: Ag34892. Range: 5-40 ng/mL.