

For Research Use Only

ATP1A1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83191-4-PBS

Basic Information

Catalog Number: 83191-4-PBS	GenBank Accession Number: BC050359	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 476	CloneNo.: 240160F9
Source: Rabbit	ENSEMBL Gene ID: ENSG00000163399	
Isotype: IgG	UNIPROT ID: P05023	
Immunogen Catalog Number: AG5763	Full Name: ATPase, Na ⁺ /K ⁺ transporting, alpha 1 polypeptide	
	Calculated MW: 113 kDa	

Applications

Tested Applications:
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00275-1, Human ATP1A1 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83191-1-PBS. Detection antibody: 83191-4-PBS. Standard: Ag5763. Range: 1.25-20 ng/mL.