

For Research Use Only

SLC6A12 Recombinant antibody, PBS Only



Catalog Number: 83195-5-PBS

Basic Information

Catalog Number: 83195-5-PBS	GenBank Accession Number: BC126215	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6539	CloneNo.: 240021C10
Source: Rabbit	UNIPROT ID: P48065	
Isotype: IgG	Full Name: solute carrier family 6 (neurotransmitter transporter, betaine/GABA), member 12	
Immunogen Catalog Number: AG18778	Calculated MW: 614 aa, 69 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83195-5-RR (SLC6A12 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83195-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83195-5-RR against Human SLC6A12 were performed. The affinity constant is below 1 pM.