

For Research Use Only

CD80/B7-1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83198-1-PBS

Basic Information

Catalog Number: 83198-1-PBS	GenBank Accession Number: NM_005191	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 941	CloneNo.: 230522D2
Source: Rabbit	UNIPROT ID: P33681	
Isotype: IgG	Full Name: CD80 molecule	

Applications

Tested Applications:
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00284-1, HUMAN CD80/B7-1 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83198-2-PBS. Detection antibody: 83198-1-PBS. Standard: Eg0060. Range: 0.625-20 ng/mL.