

For Research Use Only

CD80/B7-1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83198-2-PBS

Basic Information

Catalog Number:

83198-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_005191

GeneID (NCBI):

941

UNIPROT ID:

P33681

Full Name:

CD80 molecule

Purification Method:

Protein A purification

CloneNo.:

230522D10

Applications

Tested Applications:

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00284-1, HUMAN CD80/B7-1 Recombinant Matched Antibody Pair - BSA and Azide Free measured by Cytometric bead array. Capture antibody: 83198-2-PBS. Detection antibody: 83198-1-PBS. Standard: Eg0060. Range: 0.625-20 ng/mL.