

For Research Use Only

IGFBP4 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83200-2-PBS

Basic Information

Catalog Number:

83200-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001552.3

GeneID (NCBI):

3487

UNIPROT ID:

P22692-1

Full Name:

insulin-like growth factor binding protein 4

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

230543E12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00279-1, Human IGFBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83200-2-PBS. Detection antibody: 83200-1-PBS. Standard: Eg0825. Range: 1.56-50 ng/mL