

For Research Use Only

ATP6 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83201-1-PBS

Basic Information

Catalog Number: 83201-1-PBS	GenBank Accession Number: J01415	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4508	CloneNo.: 240041G2
Source: Rabbit	Full Name: ATP synthase 6; ATPase subunit 6	
Isotype: IgG	Calculated MW: 25 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00307-1, ATP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83201-2-PBS. Detection antibody: 83201-1-PBS. Standard: SY01389. Range: 5-80 ng/mL