

For Research Use Only

Sestrin 2 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83202-1-PBS

Basic Information

Catalog Number: 83202-1-PBS	GenBank Accession Number: BC013304	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 83667	CloneNo.: 240034F2
Source: Rabbit	UNIPROT ID: P58004	
Isotype: IgG	Full Name: sestrin 2	
Immunogen Catalog Number: AG1247	Calculated MW: 54 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00278-2, Sestrin2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83202-2-PBS. Detection antibody: 83202-1-PBS. Standard: Ag1247. Range: 1.25-20 ng/mL.