

For Research Use Only

# KLF5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83205-3-PBS

## Basic Information

<b>Catalog Number:</b> 83205-3-PBS	<b>GenBank Accession Number:</b> BC042131	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 688	<b>CloneNo.:</b> 230509D7
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q13887	
<b>Isotype:</b> IgG	<b>Full Name:</b> Kruppel-like factor 5 (intestinal)	
<b>Immunogen Catalog Number:</b> AG15267	<b>Calculated MW:</b> 457 aa, 51 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00155-1, KLF5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83205-3-PBS. Detection antibody: 83205-1-PBS. Standard: Ag15267. Range: 0.625-20 ng/mL